


MAATTABEL


JOBST® FARROWWRAP® ARMVERBAND CONFECTIE



Meetpunt (in cm)		S	M	L
Ⓒ	Omvang oksel	22-31	29-39	32-45
Ⓔ	Omvang elleboog	20-27	25-34	30-40
Ⓒ	Omvang pols	14-18	16-21	19-25
Lengte (in cm)		Kort	Normaal	Lang
I(C-G)		40-43	44-47	48-51

MAATTABEL

JOBST® FARROWWRAP® HANDVERBAND




Meetpunt (in cm)		XS	S	M	L
Ⓐ	Omvang aanzet vingers	16-18	18-20	20-22	22-24
Ⓑ	Omvang handpalm	17-19	19-21	21-23	23-26
Ⓒ	Pols	15-18	17-20	19-22	22-24

MAATTABEL


JOBST® FARROW VINGERVERBAND

Meetpunt (in cm)		1	2	3	4	5	6	7	8
Ⓐ	Omvang aanzet vingers	16,6-18,0	16,6-18,0	18,0-19,5	18,0-19,5	18,0-19,5	19,5-21,5	19,5-21,5	21,5-23,1
Ⓒ	Pols	15,6-16,7	16,7-18,1	15,9-17,0	17,0-18,5	18,5-20,4	17,2-18,7	18,7-20,6	18,3-20,1



MAATTABEL


JOBST[®] FARROWWRAP[®] TEENVERBAND



Meetpunt (in cm)		XS	S	M	L
Ⓐ	Omvang voorvoet	20-22	22-24	24-26	26-28

MAATTABEL


JOBST[®] FARROWWRAP[®] VOETVERBAND



Meetpunt (in cm)		XS	S	M	L	XL
Ⓐ	Omvang voorvoet	22-24	25-27	28-30	31-34	35-40
l Ⓐ	Voetlengte normaal	16-17	17-18	19-20	20-21	22-23
l Ⓐ	Voetlengte lang	18-19	19-20	21-22	22-23	24-25

MAATTABEL


JOBST® FARROWHYBRID® AD LINER



Meetpunt	S	M	L	XL
Schoenmaat	35-38	38-41	41-44	44-46

MAATTABEL


JOBST[®] FARROWWRAP[®] ONDERBEENVERBAND



Meetpunt (in cm)		XS	S	M	L	XL
Ⓒ	Omvang kuit	36-43	42-50	48-58	53-63	58-68
Ⓑ	Omvang enkel	21-25	25-30	30-36	36-42	42-50
I (A) - Ⓓ	Normaal	33-37	35-39	37-41	39-43	39-43
I (A) - Ⓓ	Lang	38-41	40-43	42-45	44-47	44-47

MAATTABEL


JOBST® FARROWWRAP® KNIEVERBAND



Meetpunt (in cm)		XS	S	M	L	XL	XXL
Ⓔ	Omvang knie	47-52	53-59	60-66	67-73	74-80	81-87
Ⓒ	Omvang kuit	36-43	42-50	48-58	53-63	58-68	63-73

MAATTABEL


JOBST® FARROWWRAP® BOVENBEENVERBAND



Meetpunt (in cm)		XS	S	M	L	XL
Ⓞ	Omvang bovenbeen	60-70	65-75	70-80	75-90	85-100
ⓔ	Omvang knie	47-52	53-59	60-66	67-73	74-80
Lengte (in cm)		Kort	Normaal	Lang		
I/E-G		16-20	21-26	27-31		

MAATTABEL

JOBST[®] FARROWWRAP[®] 4000



Meetpunt (in cm)		XS	S	M	L
Ⓒ	Omvang kuit	25-38	32-48	37-56	43-60
Ⓑ	Omvang enkel	17-23	19-28	25-37	29-43
I Ⓑ - Ⓓ	Normaal	30-33	30-33	30-33	30-33
I Ⓑ - Ⓓ	Lang	34-37	34-37	34-37	34-37